



## Synthesis of *cis*- and *trans*-copper(II) aminoethanoate

### Student worksheet

---



Glycine  
(aminoethanoic  
acid)



Glycinate ion  
(aminoethanoate  
ion)

Soluble copper(II) salts dissolve in water to give a solution of hexaaquacopper(II) ions,  $[\text{Cu}(\text{H}_2\text{O})_6]^{2+}$ .

Glycine (systematic name: aminoethanoic acid) has the molecular formula  $\text{C}_2\text{H}_5\text{O}_2$  and structural formula  $\text{H}_2\text{CH}_2\text{COO}$ .

When dissolved in an alkali such as sodium hydroxide, the carboxylic acid group ionises to give an ion with a single negative charge. It is called a glycinate ion (systematic name: aminoethanoate,  $\text{H}_2\text{CH}_2\text{COO}^-$ ).

The copper(II) complex of glycine can be made by adding a hot aqueous glycine solution to a hot solution of copper ethanoate in a water-ethanol mixture.

### Equipment and materials for both preparations

- Hot water bath or electric heating mantle
- 25 cm<sup>3</sup> measuring cylinder
- 2 x 100 cm<sup>3</sup> conical flasks
- 2 x 250 cm<sup>3</sup> conical flasks
- Water bath
- -10 – 110 °C thermometer
- Filter funnel and fluted filter paper
- Watch glass (2)
- 50 cm<sup>3</sup> round-bottom flask (Quickfit)
- Condenser (Quickfit) to fit the 50 cm<sup>3</sup> round-bottom flask
- Bunsen Burner, tripod and gauze
- Sample tube (2)
- Ethanol
- Copper(II) ethanoate-1-water
- Glycine
- Ice

### Synthesis of *cis*- $\text{Cu}(\text{H}_2\text{CH}_2\text{COO})_2 \cdot \text{H}_2\text{O}$

#### Method

**Care:** Wear eye protection. Copper(II) ethanoate-1-water is harmful and an irritant. Ethanol is flammable. Work away from naked flames. As industrial denatured alcohol or ethanol formulated with methanol (IDA/IMS – industrial methylated spirits) is being used, there may be toxicity issues from methanol (or other adulterants) so work in a well-ventilated lab.

1. Using a measuring cylinder, measure out 25 cm<sup>3</sup> of ethanol into a 100 cm<sup>3</sup> conical flask and place the flask in a hot water bath at 65 °C.
2. Using a measuring cylinder, measure 25 cm<sup>3</sup> of water into a 250 cm<sup>3</sup> conical flask. Weigh out 2.0 g (0.01 mol) of copper(II) ethanoate-1-water and add to the flask. Place the flask in the water bath and stir the mixture until the solid dissolves. When the solution reaches 65 °C pour the warm ethanol into the copper(II) ethanoate solution.
3. Using a measuring cylinder, measure 25 cm<sup>3</sup> of water into a 100 cm<sup>3</sup> conical flask. Weigh out 1.5 g (0.02 mol) of glycine and add to the flask. Place the flask in the water bath. Stir the mixture until the solid dissolves. When the solution reaches 65 °C pour it into the flask containing the copper(II) ethanoate solution.

- Cool the mixture in an ice bath. Light blue needles should crystallise. Filter the cold mixture through a fluted filter paper. Collect the filtrate in clean 250 cm<sup>3</sup> conical flask and keep it for the synthesis of *trans*-Cu(H<sub>2</sub>NCH<sub>2</sub>COO)<sub>2</sub>·H<sub>2</sub>O.
- Open up the filter paper and lay it on a watch glass. Cover the crystals with a piece of clean filter paper and leave it to dry in a fume cupboard at room temperature.
- Label a sample tube with the name of the product, your name and the date. Weigh the labelled sample tube and record its mass.
- Scrape the dried crystals into the weighed sample tube. Weigh the tube again. Record its mass.

### Calculations

Calculate the theoretical yield and the percentage yield.

### Synthesis of *trans*-Cu(H<sub>2</sub>CH<sub>2</sub>COO)<sub>2</sub>·H<sub>2</sub>O

#### Method

- Weigh about 1.5 g of *cis*-Cu(H<sub>2</sub>NCH<sub>2</sub>COO)<sub>2</sub>·H<sub>2</sub>O and 1 g of glycine into a 50 cm<sup>3</sup> round-bottom flask.
- Use a measuring cylinder to pour 10 cm<sup>3</sup> of the filtrate from the synthesis of *cis*-Cu(H<sub>2</sub>NCH<sub>2</sub>COO)<sub>2</sub>·H<sub>2</sub>O into the flask, fit a condenser and reflux the mixture for one hour (use a hot water bath or electric heating mantle).
- Filter the hot mixture through a fluted filter paper to obtain blue-violet platelets of *trans*-Cu(H<sub>2</sub>NCH<sub>2</sub>COO)<sub>2</sub>·H<sub>2</sub>O.

Note that if the mixture all dissolves during the reflux, then the *cis*-product will re-precipitate not the *trans* isomer.

- Open up the filter paper and lay it on a watch glass. Cover the crystals with a piece of clean filter paper and leave it to dry in a fume cupboard at room temperature.
- Label a sample tube with the name of the product, your name and the date.
- Scrape the dried crystals into the weighed sample tube.

It is not possible to get a meaningful value for the percentage yield as the exact quantity of *cis*-Cu(H<sub>2</sub>NCH<sub>2</sub>COO)<sub>2</sub>·H<sub>2</sub>O is not known.

### Finally

Why does the structure change from the *cis* isomer to the *trans* isomer? The initial reaction mixture contains both *cis* and *trans* isomers in equilibrium, but the *cis* isomer crystallises faster. However the *trans* isomer is thermodynamically more stable and less soluble, so on continued heating, the only product is the *trans* isomer.

- Make molecular models and draw structures for the two geometrical isomers of copper(I) aminoethanoate.
- Do you think there are any benefits of using one isomer rather than the other as a source of the micronutrient copper?