

# **Preparation of zinc sulfate**

#### Student worksheet

## Making zinc sulfate

Zinc sulfate is a complex fertiliser. This simple salt is a source of the micronutrient zinc and the secondary nutrient sulfur. However, it is as a source of zinc that it is important.

Zinc sulfate is produced by reacting zinc carbonate with sulfuric acid:

 $ZnCO_3(s) + H_2SO_4(aq) \rightarrow ZnSO_4(aq) + H_2O(l) + CO_2(g)$ 

Zinc sulfate used for pharmaceutical preparations is made from the reaction of sulfuric acid with high purity zinc oxide:

 $ZnO(s) + H_2SO_4(aq) \rightarrow ZnSO_4(aq) + H_2O(l)$ 

Zinc sulfate crystallises from aqueous solution as a heptahydrate, zinc sulfate-7-water,  $ZnSO_4.7H_2O$ .

You can make zinc sulfate-7-water in the laboratory by reacting zinc carbonate with dilute sulfuric acid.

## **Equipment and materials**

- Weighing bottle (or small beaker)
- 250 cm<sup>3</sup> beaker
- Bunsen burner, tripod and gauze
- 25 cm<sup>3</sup> measuring cylinder
- Hot water bath
- Filter funnel and filter paper

- Thermometer  $(0 110 \degree C)$
- Evaporating basin
- Spatula
- Sample bottle
- 1 mol dm<sup>-3</sup> sulfuric acid
- Small granules of zinc carbonate

# Method

**Care:** Wear eye protection. 1 mol dm<sup>-3</sup> sulfuric acid is an irritant.

- 1. Weigh out about 3.5 g of zinc carbonate into a weighing bottle or small beaker.
- 2. Using a measuring cylinder, measure out 25 cm<sup>3</sup> of 1.0 mol dm<sup>-3</sup> sulfuric acid into a 250 cm<sup>3</sup> beaker.
- 3. Use a spatula to add zinc carbonate to the acid. Add the solid a little at a time, stirring between additions and allowing any effervescence to die away before adding more.
- 4. Filter the solution into an evaporating basin. This removes any unreacted zinc carbonate from the mixture.
- 5. Put the evaporating basin on a tripod and gauze. Slowly evaporate the solution until it is about one-fifth of its original volume. Caution: Do not boil the solution as it may spit.
- 6. Allow the concentrated solution to cool until crystals form.
- 7. Filter off the crystals and put the filter paper and crystals on a watch glass and dab dry with another piece of filter paper. Cover them with a piece of clean filter paper and leave them to dry at room temperature.





- 8. Label a sample tube with the name of the product, your name and the date. Weigh the labelled sample tube and record its mass.
- 9. Tip your dry product into the sample tube. Weigh the tube again. Record its mass.

#### Calculations

You used an excess of zinc carbonate and so the theoretical yield depends on the volume of 1 mol  $dm^{-3}$  sulfuric acid used.

Calculate

- 1. the theoretical yield of zinc sulfate-7-water;
- 2. the percentage yield of zinc sulfate-7-water.