

Iron by thiocyanate assay

Teacher and technician worksheet

Equipment and materials

Each student or pair of students will require:

- burettes x 3 (for greater accuracy, graduated pipettes could be used)
- 10 cm³ pipette
- 100 cm³ beaker x 7
- colorimeter and suitable filter (blue) A solution of the complex displays maximum absorption at 480 nm
- iron(III) ammonium sulfate solution containing 0.050 g dm⁻³ Fe³⁺ (50 ppm) (15 cm³)
- ammonium thiocyanate solution, 1 mol dm⁻³ (70 cm³)
- solution of unknown Fe³⁺ concentration (10 cm³)

Make sure students wear eye protection.

Solution preparations

- Iron(III) ammonium sulfate solution, 0.050 gdm⁻³ Fe³⁺ (50 ppm)
Measure 0.432 g of finely powdered iron(III) ammonium sulfate-12-water into a 100 cm³ beaker, add 200 cm³ of 2 mol dm⁻³ sulfuric(VI) acid. Set this on a stirrer to dissolve. You may have to leave this over night to fully dissolve, but the stirrer does not need to be on. Add the contents and the washings to a 1 dm³ volumetric flask. Make up to volume with distilled or deionised water.
- Ammonium thiocyanate solution, 1 mol dm⁻³
Weigh 38 g of solid ammonium thiocyanate or 49 g of potassium thiocyanate into a 500 cm³ beaker. Add about 200 cm³ of distilled or deionised water and stir to dissolve the solid. Transfer the solution to a 500 cm³ volumetric flask and make up to the mark with distilled or deionised water.