

## Determination of thiocyanate using iron(III)

### Teacher and technician worksheet

---

#### Equipment and materials

Each student or pair of students will require:

- burettes x 3
- 100 cm<sup>3</sup> volumetric flasks x 7 (or use one, thoroughly washing it between samples)
- colorimeter and suitable filter (blue) - a solution of the complex displays maximum absorption at 480 nm
- potassium thiocyanate containing 250 mg dm<sup>-3</sup> thiocyanate (250 ppm) (30 cm<sup>3</sup>)
- iron(III) chloride solution 0.41 mol dm<sup>-3</sup> (70 cm<sup>3</sup>)
- solution of unknown thiocyanate concentration (10 cm<sup>3</sup>)

Make sure students wear eye protection. Iron(III) chloride solution is an irritant.

#### Solution preparations

- Standard solution of potassium thiocyanate containing 250 mg dm<sup>-3</sup>  
Dissolve 4.5 g potassium thiocyanate in deionised water and dilute to 500 cm<sup>3</sup> in a volumetric flask. This solution contains 5000 mg dm<sup>-3</sup> thiocyanate ions. Dilute 50 cm<sup>3</sup> of this solution to 1 dm<sup>3</sup> in a volumetric flask. The diluted solution contains 250 mg dm<sup>-3</sup> (250 ppm) thiocyanate ions.
- Iron(III) chloride solution, 0.41 mol dm<sup>-3</sup>  
Dissolve 50 g iron(III) chloride-6-water in a beaker in about 200 cm<sup>3</sup> 1 mol dm<sup>-3</sup> hydrochloric acid. Make up to 250 cm<sup>3</sup> in a volumetric flask with more 1 mol dm<sup>-3</sup> hydrochloric acid.