

Determining the pK_a of 2-hydroxybenzoic acid

Student worksheet

Health and safety note

Wear eye protection and ensure no naked flames. 0.10 mol dm^{-3} sodium hydroxide solution is an irritant. 95% ethanol is highly flammable. 2-hydroxybenzoic acid and 95% ethanol are both harmful.

Principle

2-hydroxybenzoic acid is a weak acid






Figure 1 2-hydroxybenzoic acid ionises in aqueous solution. It is a weak acid.

This may be represented:



or, more simply



Its acid dissociation constant, K_a , is given by:

$$K_a = \frac{[\text{H}_3\text{O}^+][\text{A}^-]}{[\text{HA}]}$$

often written simply as:

$$K_a = \frac{[\text{H}^+][\text{A}^-]}{[\text{HA}]}$$

The pK_a value is given by:

$$pK_a = -\log_{10} K_a$$

Taking logarithms, the following relationships are derived

$$\text{pH} = pK_a - \log \frac{[\text{HA}]}{[\text{A}^-]}$$

or

$$\text{pH} = pK_a + \log \frac{[\text{A}^-]}{[\text{HA}]}$$

When the acid is 'half-neutralised',

$[\text{A}^-] = [\text{HA}]$, and $\text{pH} = pK_a$.





Figure 2 A graph of pH for the titration of a weak acid with sodium hydroxide solution.

By titrating a solution of the 2-hydroxybenzoic acid against a strong alkali, such as sodium hydroxide solution, the pH at the half-way point can be determined and this gives the pK_a of aspirin.

2-hydroxybenzoic acid and sodium hydroxide react in a 1:1 mole ratio:



Equipment and materials

- Balance
- 50 cm³ burette
- 250 cm³ beaker
- Glass stirring rod
- 10 cm³ and 100 cm³ measuring cylinders
- Spatula
- pH probe and pH meter
- 2-hydroxybenzoic acid – Harmful
- 95% ethanol – Highly flammable, Harmful
- 0.10 mol dm⁻³ sodium hydroxide solution – Irritant

Method

1. Fill a burette with 0.10 mol dm⁻³ sodium hydroxide solution.
2. Weigh 0.28 g of 2-hydroxybenzoic acid into a 250 cm³ beaker. Add 10 cm³ of 95% ethanol, Stir with a glass rod and when the solid has dissolved add 90 cm³ of deionised water. Stir the mixture until it is homogeneous.
3. Place a pH probe in the solution and connect it to a pH meter.
Note: The pH probe should have been calibrated using suitable buffer solutions.
4. Add 2 cm³ quantities of sodium hydroxide solution from the burette to the beaker, stirring well between additions and recording the pH.
5. Near the end-point the pH begins to rise rapidly. So after you have added 18 cm³ of sodium hydroxide solution begin adding it in 0.5 cm³ portions. After about 22 cm³ start adding in 2 cm³ portions again. Continue until total of 36 cm³ has been added.

Processing data

1. Plot a graph of pH against volume of 0.10 mol dm⁻³ sodium hydroxide solution added.
2. From the graph, calculate the end-point of the titration.
3. Check this against the expected value by calculating the number of moles of 2-hydroxybenzoic acid used (relative molecular mass of 2-hydroxybenzoic acid = 138) and, therefore, the volume of 0.10 mol dm⁻³ sodium hydroxide solution needed to react with it in a 1:1 mole ratio.
4. At the half-way point to the end point, $[\text{HOC}_6\text{H}_4\text{COO}^-] = [\text{H}^+]$.
From the graph, estimate the pH at the half-way point of the titration and, therefore, a value for the pK_a of 2-hydroxybenzoic acid.
5. Calculate K_a of 2-hydroxybenzoic acid.