

## Iron by 1,10-phenanthroline assay

### Teacher and technician worksheet

---

#### Equipment and materials

Each student or pair of students will require:

- burettes x 3
- 1 cm<sup>3</sup> pipette
- 2 cm<sup>3</sup> pipette
- 100 cm<sup>3</sup> volumetric flask x 7 (or re-use one flask)
- colorimeter and suitable filter (green) – a solution of the complex displays maximum absorption at 510 nm
- iron(II) ammonium sulfate solution containing 0.100 g dm<sup>-3</sup> Fe<sup>2+</sup> (100 ppm), (15 cm<sup>3</sup>)
- hydroxylamine hydrochloride solution (if used), 1.5 mol dm<sup>-3</sup> (7 cm<sup>3</sup>)
- 1,10-phenanthroline solution, 5 x 10<sup>-3</sup> mol dm<sup>-3</sup> (105 cm<sup>3</sup>)
- sodium ethanoate solution, 1 mol dm<sup>-3</sup> (14 cm<sup>3</sup>)
- solution of unknown Fe<sup>2+</sup> concentration, (15 cm<sup>3</sup>)

Make sure students wear eye protection. Hydroxylamine hydrochloride solution is harmful.

#### Solution preparation

- Iron(II) ammonium sulfate solution, 0.100 g dm<sup>-3</sup> Fe<sup>2+</sup> (100 ppm)  
Weigh out 0.702 g of iron(II) ammonium sulfate-6-water, Fe(NH<sub>4</sub>)<sub>2</sub>(SO<sub>4</sub>)<sub>2</sub>·6H<sub>2</sub>O, dissolve in deionised water and make up to 1 dm<sup>3</sup>.
- 1,10-phenanthroline solution, 5 x 10<sup>-3</sup> mol dm<sup>-3</sup>  
Weigh out 0.90 g of 1,10-phenanthroline, dissolve in deionised water, warming if necessary but do not allow the solution to boil. Discard the solution if it darkens. Make up to 1 dm<sup>3</sup>.
- Sodium ethanoate solution, 1 mol dm<sup>-3</sup>  
Weigh out 13.6 g of sodium ethanoate-3-water, dissolve in deionised water and make up to 100 cm<sup>3</sup>.
- Hydroxylamine hydrochloride solution, 1.5 mol dm<sup>-3</sup>  
Weigh out 10.4 g of hydroxylamine hydrochloride, dissolve in deionised water and make up to 100 cm<sup>3</sup>.