


Iron by 1,10-phenanthroline assay

Student worksheet


Principle

Iron(II) ions form a deep red solution with 1,10-phenanthroline. Generally, it is suitable for solutions that are 5 mg dm^{-3} , or less, $\text{Fe}^{2+}(\text{aq})$. The colour is due to the complex that forms:



You can use this reaction for the quantitative analysis of low concentrations of $\text{Fe}^{2+}(\text{aq})$ in solution. You can find the concentration of the solution of Fe^{2+} using a colorimeter. You can also use simple colour matching although the results will be less precise. This method is suitable for solutions that contain 5 mg dm^{-3} , or less, $\text{Fe}^{2+}(\text{aq})$.

If iron is present in the +3 oxidation state, it must be reduced to the +2 oxidation state before complex formation. This can be done by reaction with hydroxylamine hydrochloride:



Equipment and materials

- burettes x 3
- 1 cm^3 pipette
- 2 cm^3 pipette
- 100 cm^3 volumetric flask x 7 (or re-use one flask)
- colorimeter and suitable filter (green) – a solution of the complex displays maximum absorption at 510 nm
- iron(II) ammonium sulfate solution containing $0.100 \text{ g dm}^{-3} \text{ Fe}^{2+}$ (100 ppm), (15 cm^3)
- hydroxylamine hydrochloride solution (if used), 1.5 mol dm^{-3} (7 cm^3)
- 1,10-phenanthroline solution, $5 \times 10^{-3} \text{ mol dm}^{-3}$ (105 cm^3)
- sodium ethanoate solution, 1 mol dm^{-3} (14 cm^3)
- solution of unknown Fe^{2+} concentration, (15 cm^3)

Method

Care: Wear eye protection. Hydroxylamine hydrochloride solution is harmful.

1. Fill two burettes, one with the standard iron(II) solution and one with deionised water.
2. Label six 50 or 100 cm^3 beakers A to F and use the burettes to add the solutions shown in the table:

Volumetric flask	A	B	C	D	E	F
Volume of standard iron(II) solution / cm ³	5.0	4.0	3.0	2.0	1.0	0.0

3. To each volumetric flask, use a 1 or 2 cm³ pipette to add 1.0 cm³ of 1.5 mol dm⁻³ hydroxylamine hydrochloride solution (if used) and 2.0 cm³ of 1 mol dm⁻³ sodium ethanoate solution. From a burette add 15.0 cm³ of 5 x 10⁻³ mol dm⁻³ 1,10-phenanthroline solution. Make up the volume to 100 cm³ using distilled or deionised water.

Beaker	A	B	C	D	E	F
ppm iron	5	4	3	2	1	0

4. Mix each solution well and leave to stand for 10 minutes.
5. Measure the absorbance of each solution. If the colour of the solutions is too intense for the colorimeter to measure, dilute the solution by adding more water to the beakers from the burette and measure the absorbance of the diluted solutions.
6. Plot a graph of absorbance (y axis) against Fe²⁺(aq) concentration (in ppm iron) (x axis) for beakers A-F.
7. Using a burette, add 10 cm³ of the iron(II) solution of unknown Fe²⁺ concentration to a 100 cm³ volumetric flask. Use a 1 or 2 cm³ pipette to add 1.0 cm³ of 1.5 mol dm⁻³ hydroxylamine hydrochloride solution (if used) and 2.0 cm³ of 1 mol dm⁻³ sodium ethanoate solution. From a burette add 15.0 cm³ of 5 x 10⁻³ mol dm⁻³ 1,10-phenanthroline solution. Make up the volume to 100 cm³ using distilled or deionised water.
8. Measure the absorbance of the solution made up using of the iron(II) solution whose concentration is unknown.
9. Use the graph to find the concentration of Fe²⁺(aq) as ppm iron in the unknown. Remember that the iron(II) solution has been diluted ten-fold in the volumetric flask.