
Some reactions of nitrogen dioxide - student sheet

Procedure

- 1. Cover the worksheet with a clear plastic sheet.
- 2. Place the base of the Petri dish directly over the circle below. Place the reaction vessel in the centre of the Petri dish.
- 3. Put two drops of full-range indicator where shown.
- 4. At another corner of the triangle, place two drops of ammonia solution. Place the lid on the Petri dish and wait for the indicator drop to change colour.
- 5. Remove the lid from the Petri dish and, using a piece of tissue, mop up the drop of ammonia.
- 6. At the two remaining corners of the triangle, add the two other test solutions.
- 7. Add a few copper metal turnings to the reaction vessel, followed by three drops of nitric acid. Quickly replace the lid on the Petri dish.
- 8. Record all your observations over the next 15 min.





This resource was downloaded from https://rsc.li/3dUhlKN

Questions

1. What explanations can you give for your observations?

Health, safety and technical notes

- Read our standard health and safety guidance here https://rsc.li/3LVOUbN
- Students must wear suitable eye protection (Splash resistant goggles to BS EN166 3).
- Nitrogen dioxide is extremely toxic and corrosive if inhaled, with sometimes delayed effects. It is important to ensure that the amount of NO₂ generated does not result in significant leakage from the Petri dish (see CLEAPSS Hazcard <u>HC068b</u>). No more than 2 copper turnings should be used.
- Concentrated Nitric acid, HNO₃(aq), 5 mol dm⁻³, is CORROSIVE and gives of toxic fumes (see CLEAPSS Hazcard <u>HC067</u>).
- Potassium iodate(V), KIO₃(aq),0.1 mol dm⁻³, Potassium bromate(V), KBrO₃(aq), 0.1 mol dm⁻³ (see CLEAPSS Hazcard <u>HC080</u>), Potassium bromide, KBr(aq), 0.2 mol dm⁻³ and Potassium iodide, KI(aq), 0.2 mol dm⁻³ are low hazard (see CLEAPSS Hazcard <u>HC047b</u>).
- Ammonia solution, NH₃(aq), 3 mol dm⁻³ is corrosive and a respiratory IRRITANT (see CLEAPSS Hazcard <u>HC006</u>).

