


## Measuring the amount of vitamin C in fruit drinks – student sheet


In this experiment you will be finding out how much vitamin C there is in a fruit drink.

The chemical name for vitamin C is ascorbic acid. The basis of the experiment is as follows.


A known amount of iodine is generated by the reaction between iodate, iodide and sulfuric acid:



A measured amount of fruit drink is added. The ascorbic acid in the drink reacts quantitatively with some of the iodine as the iodine is in excess:



The excess iodine is then titrated against standard sodium thiosulfate solution:



From the titration results the amount of iodine that reacts with the sodium thiosulfate solution can be found.

Since the total amount of iodine originally formed is known the amount that reacts with the ascorbic acid is found by difference.

Therefore, the amount of ascorbic acid that reacts with this amount of iodine can be found.

### Procedure

1. Set up the microscale titration apparatus.
2. Fill the apparatus with sodium thiosulfate solution.
3. Using the glass pipette, add 2 cm<sup>3</sup> of potassium iodate solution to the beaker.
4. Measure, using the measuring cylinder, 3 cm<sup>3</sup> of potassium iodide solution, then add this to the beaker. (Note: the potassium iodide solution is added in slight excess.)
5. Add three drops of sulfuric acid. A yellow-brown colour appears due to iodine.
6. Add a few drops of starch solution. A deep blue-black colour forms.
7. Using the glass pipette, add 1 cm<sup>3</sup> of the fruit drink to the beaker and swirl gently.
8. Titrate the remaining iodine in the beaker against the sodium thiosulfate solution. (The beaker can be swirled very gently to mix the chemicals. Alternatively, the tip of a plastic pipette can be used as a mini stirring rod.) The disappearance of the deep blue-black colour marks the end-point.
9. Do a duplicate titration and check the agreement between the two titres. If it is acceptable, take the mean value of the two titres and use it for your calculations.

## Calculations

A specimen result and calculation is given below.

Study this carefully and use it as a guide for working out the vitamin C content of your fruit drink.

The volume of thiosulfate delivered during the titration =  $0.74 \text{ cm}^3$ .

The concentration of thiosulfate =  $0.010 \text{ mol dm}^{-3}$ .

Therefore, the number of moles of thiosulfate =

$$\frac{0.74 \times 0.010}{1000} = 7.4 \times 10^{-6}$$

Therefore, the number of moles of iodine that this reacts with during the titration is  $3.7 \times 10^{-6}$

The total number of moles of iodine produced in the reaction between iodate, iodide and sulfuric acid based on using  $2 \text{ cm}^3$  of iodate with a concentration of  $0.0012 \text{ mol dm}^{-3}$  =

$$\frac{3 \times 2 \times 0.0012}{1000} = 7.2 \times 10^{-6}$$

Therefore, the number of moles of iodine that reacts with the ascorbic acid is  $7.2 \times 10^{-6} - 3.7 \times 10^{-6} = 3.5 \times 10^{-6}$ .

Since 1 mole of iodine reacts with 1 mole of ascorbic acid then the number of moles of ascorbic acid is also  $3.5 \times 10^{-6}$ . The volume of the fruit juice used is  $1 \text{ cm}^3$ .

Therefore, the number of moles of ascorbic acid in  $1000 \text{ cm}^3 = 3.5 \times 10^{-3}$ . The relative molar mass of ascorbic acid = 174.12 g.

Therefore, the mass of ascorbic acid (in  $1000 \text{ cm}^3$ ) =  $174.12 \times 3.5 \times 10^{-3} = 0.609 \text{ g}$ .

Therefore, the vitamin C content of the fruit drink = 61mg per  $100 \text{ cm}^3$ .

## Health, safety and technical notes

- Read our standard health and safety guidance here <https://rsc.li/3SRKCow>
- Wear eye protection.
- Sulfuric acid,  $1 \text{ mol dm}^{-3}$  is a skin/eye irritant (see CLEAPSS Hazcard [HC098a](#)).
- Sodium thiosulfate,  $0.010 \text{ mol dm}^{-3}$ , potassium iodate,  $0.001 \text{ mol dm}^{-3}$  and potassium iodide,  $0.005 \text{ mol dm}^{-3}$  solutions are of low hazard, as are the starch solution and fruit juices (see CLEAPSS Hazcards [HC095a](#), [HC080](#), [HC047b](#)).