

Making nylon - the 'nylon rope trick'

Introduction

A solution of decanedioyl dichloride in cyclohexane is floated on an aqueous solution of 1,6-diaminohexane. Nylon forms at the interface and can be pulled out as fast as it is produced forming a long thread – the 'nylon rope'.

Equipment

Apparatus

- Eye protection
- Beaker, 25 cm³
- A pair of tweezers
- Retort stand
- Boss and clamp

Chemicals

- 1,6-diaminohexane (hexamethylene diamine, hexane-1,6-diamine, H₂N(CH₂)₆(NH₂), 2.2 g
- Decanedioyl dichloride (sebacoyl chloride, ClOC(CH₂)₈COCl), 1.5 g
- Cyclohexane, 50 cm³
- Deionised water, 50 cm³

Health, safety and technical notes

- Read our standard health and safety guidance here <https://rsc.li/3TnTkun>
- Always wear eye protection.
- 1,6-diaminohexane (hexamethylene diamine, hexane-1,6-diamine, H₂N(CH₂)₆(NH₂) is corrosive, harmful if swallowed or in contact with skin and a respiratory irritant (see CLEAPSS Hazcard [HC003b](#)).
- Decanedioyl dichloride (sebacoyl chloride, ClOC(CH₂)₈COCl) is corrosive and harmful if swallowed (see CLEAPSS Hazcard [HC041](#)).
- Cyclohexane is highly flammable, and a skin/respiratory irritant (see CLEAPSS Hazcard [HC045b](#)).

Disposal

Dispose of the mixture as follows:

- First shake the reaction to mix the two layers. A lump of nylon will be produced which can be removed with tweezers, rinsed well with water, and disposed as solid waste.
- The remaining liquids can be mixed with detergent and washed down the sink.
- Failure to do this may result in the polymerisation taking place in the sink, leading to a blockage.
- The organic liquid can be disposed of in a hydrocarbon waste bottle stored in a flammable cabinet.

Procedure

Work in a well-ventilated laboratory. Wear eye protection and disposable nitrile gloves when pulling out the thread.


Before the demonstration

1. Make up a solution of 2.2 g of 1,6-diaminohexane in 50 cm³ of deionised water. This solution is approximately 0.4 mol dm⁻³.

2. Make up a solution of 1.5 g of decanedioyl dichloride in 50 cm³ of cyclohexane. This solution is approximately 0.15 mol dm⁻³.

The demonstration

1. Pour 5 cm³ of the aqueous diamine solution into a 25 cm³ beaker. Carefully pour 5 cm³ of the cyclohexane solution of the acid chloride on top of the first solution so that mixing is minimised. Do this by pouring the second solution down the wall of the beaker, or pour it down a glass rod.
2. The cyclohexane will float on top of the water without mixing.
3. Place the beaker below a stand and clamp as shown (see below). A greyish film of nylon will form at the interface.
4. Pick up a little of this with a pair of tweezers and lift it slowly and gently from the beaker. It should draw up behind it a thread of nylon.
5. Pull this over the rod of the clamp so that this acts as a pulley.
6. Continue pulling the nylon thread at a rate of about half a metre per second. It should be possible to pull out several metres.
7. Take care, the thread will be coated with unreacted monomer and may in fact be a narrow, hollow tube filled with monomer solution. Wearing disposable gloves is essential.



Notes

- The beaker is rather small so allow the audience as close as possible consistent with comfort and safety.
- Point out that this demonstration is different from the industrial method of making nylon, which takes place at a higher temperature. Molten nylon is then forced through multi-holed 'spinnerets' to form the fibres.
- The reaction is a condensation polymerisation:
- $n\text{H}_2\text{N}(\text{CH}_2)_6\text{NH}_2 + n\text{ClOC}(\text{CH}_2)_8\text{COCl} \rightarrow \text{H}_2\text{N}[(\text{CH}_2)_6\text{NHCO}(\text{CH}_2)_8]_n\text{COCl} + n\text{HCl}$
- The nylon formed is nylon 6-10 so called because of the lengths of the carbon chains of the monomers.

- Nylon 6 – 6 can be made using hexanedioyl dichloride (adipoyl chloride). The diamine is present in excess to react with the hydrogen chloride that is eliminated. An alternative procedure is to use the stoichiometric quantity of diamine dissolved in excess sodium hydroxide solution.