



A practical investigation of *cis-trans* isomerism in transition metal complexes

This resource accompanies the article **Teaching isomerism at post-16** in *Education in Chemistry* where you can find more ideas and support for this topic: rsc.li/3PvYVPw.

It is part of our **Challenging plants: fertilisers** collection of resources and experiments.

Resource components

	
<p>Student sheet: includes full instructions for the practical activities.</p>	<p>Presentation: slides include pre-lab questions and post-lab synoptic questions, with model answers.</p>

Learning objectives

- 1 Recognise why transition metal complexes can show *cis-trans* isomerism.
- 2 Identify further examples of structural and stereoisomerism in transition metal complex ions and other species.
- 3 Apply your knowledge and understanding of transition metal chemistry to the synthesis of the *cis-trans* isomers of copper(II) aminoethanoate.

In this resource, learners investigate isomerism in transition metal complexes and apply their knowledge and understanding to the synthesis of the *cis-* and *trans-*isomers of copper(II) aminoethanoate. Encourage learners to make links with other topics such as quantitative chemistry, amino acids and isomerism in organic compounds with the pre- and post-lab questions in the PowerPoint.

How to use this resource

This resource centres around the synthesis of the *cis-* and *trans-*isomers of copper(II) aminoethanoate. Owing to the time constraints of a typical lesson you can choose to complete the activity in two parts. You can store the *cis*-copper(II) aminoethanoate and the filtrate from Part 1 between lessons ready for isomerism to *trans*-copper(II) aminoethanoate in Part 2.

Prior to completing the synthesis, ask learners to complete the pre-lab questions on the linked PowerPoint. These questions are designed to help learners understand the

chemistry involved in the synthesis. Run through the answers immediately before starting the practical.

Learners synthesise copper(II) aminoethanoate by the reaction of aqueous copper(II) ions and glycine. The product contains a mixture of the *cis*- and *trans*-isomers in equilibrium. As the *cis*-isomer crystallises faster it is isolated first.

Once the *cis*-isomer is isolated, it is isomerised to produce a sample of the *trans*-isomer. The *trans*-isomer is more thermodynamically stable and is less soluble and so can be isolated from the equilibrium mixture of the two isomers after sustained heating.

The PowerPoint provides details of the synthesis and isomerism, along with scaffolding to support learners with the steps needed to calculate the percentage yield of the *cis*-isomer.

Four post-lab synoptic questions and answers are provided on the PowerPoint to encourage learners to retrieve prior knowledge and make links between the topics. Set these as homework or ask learners to complete them during Part 2 of the practical activity, while the solution is being heated under reflux.

Use the pre- and post-lab questions as a standalone activity if time or resource constraints prevent completion of the synthesis. Alternatively, get learners to synthesise the *cis*-isomer only and ask them to compare its appearance with a pre-made sample of the *trans*-isomer.

The extension activity is an open question connected to the use of copper(II) aminoethanoate as a fertiliser. Use this as a discussion as learners complete the practical or as a follow-up question. Additionally, you can ask learners to research and present the work of Edith Ellen Humphrey who synthesised isomers of octahedral complexes and was a pioneering woman in chemistry (see rsc.li/4gkt5kK).

Technical notes

Read our standard health and safety guidance (rsc.li/3zyJLkx) and carry out a risk assessment before running any live practical.



Equipment and materials for both preparations (per group)

- Splash-proof goggles
- Hot water bath or electric heating mantle
- 25 cm³ measuring cylinder
- 2 × 100 cm³ conical flasks
- 2 × 250 cm³ conical flasks
- Water bath
- -10–110°C thermometer
- Filter funnel and fluted filter paper

- 2 × watch glasses
- 50 cm³ round-bottom flask (Quickfit)
- Condenser (Quickfit) to fit the 50 cm³ round-bottom flask
- Bunsen burner, tripod and gauze
- 2 × sample tubes
- Ethanol (DANGER: flammable), 25 cm³
- Copper(II) ethanoate monohydrate (DANGER: harmful; irritant), 2.0 g
- Glycine, 2.5 g
- Distilled water, 50 cm³
- Ice

Safety and hazards

- Wear splash-proof goggles.
- Work in a well-ventilated laboratory.
- Ethanol is flammable: work away from naked flames.

Chemicals supplied for the practical	Hazards
Copper(II) ethanoate monohydrate  DANGER	Harmful if swallowed. Causes skin irritation. Causes serious eye damage. May cause respiratory irritation. Very toxic to aquatic life with long lasting effects. See CLEAPSS Hazcard HC027B .
Ethanol  DANGER	Highly flammable liquid and vapour. IDA is also: harmful if swallowed and may cause damage to organs. See CLEAPSS Hazcard HC040A .

Part 1 – synthesis of $\text{cis-Cu}(\text{H}_2\text{NCH}_2\text{COO})_2 \cdot \text{H}_2\text{O}$

Method

1. Using a measuring cylinder, measure out 25 cm³ of ethanol into a 100 cm³ conical flask and place the flask in a hot water bath at 65°C.
2. Using a measuring cylinder, measure 25 cm³ of distilled water into a 250 cm³ conical flask. Weigh out 2.0 g (0.01 mol) of copper(II) ethanoate monohydrate

and add it to the flask. Place the flask in the water bath and stir the mixture until the solid dissolves. When the solution reaches 65°C, pour the warm ethanol into the copper(II) ethanoate solution.

- Using a measuring cylinder, measure 25 cm³ of distilled water into a 100 cm³ conical flask. Weigh out 1.5 g (0.02 mol) of glycine and add it to the flask. Place the flask in the water bath. Stir the mixture until the solid dissolves. When the solution reaches 65°C, pour it into the flask containing the copper(II) ethanoate solution.
- Cool the mixture in an ice bath. Light blue needles should crystallise. Filter the cold mixture through a fluted filter paper. Collect the filtrate in a clean 250 cm³ conical flask and keep it for the synthesis of *trans*-Cu(H₂NCH₂COO)₂·H₂O.
- Open up the filter paper and lay it on a watch glass. Cover the crystals with a piece of clean filter paper and leave to dry in a fume cupboard at room temperature.
- Label a sample tube with the name of the product, your name and the date. Weigh the labelled sample tube and record its mass.
- Scrape the dried crystals into the weighed sample tube. Weigh the tube again. Record its mass.

Calculations

Ask learners to calculate the theoretical yield and the percentage yield. Scaffolding for the calculation is on Slide 9 of the PowerPoint if learners need support.

Part 2 – synthesis of *trans*-Cu(H₂NCH₂COO)₂·H₂O

Method

- Weigh 1.5 g of *cis*-Cu(H₂NCH₂COO)₂·H₂O and 1 g of glycine into a 50 cm³ round-bottom flask.
- Use a measuring cylinder to pour 10 cm³ of the filtrate from the synthesis of *cis*-Cu(H₂NCH₂COO)₂·H₂O into the flask, fit a condenser and reflux the mixture for one hour (use a hot water bath or electric heating mantle).
- Filter the hot mixture through a fluted filter paper to obtain blue-violet platelets of *trans*-Cu(H₂NCH₂COO)₂·H₂O. Note that if the mixture all dissolves during the reflux, then the *cis*-product will re-precipitate not the *trans*-isomer.
- Open up the filter paper and lay it on a watch glass. Cover the crystals with a piece of clean filter paper and leave to dry in a fume cupboard at room temperature.
- Label a sample tube with the name of the product, your name and the date.
- Scrape the dried crystals into the weighed sample tube. It is not possible to get a meaningful value for the percentage yield as the exact quantity of *cis*-Cu(H₂NCH₂COO)₂·H₂O is not known.