

Microscale universal indicator

These technician notes are part of a collection of microscale chemistry resources at: rsc.li/4iiljbl. Integrated instructions for this experiment are also available from rsc.li/4i8Muq6.

Introduction

Use this activity to introduce learners aged 11–16 to the concept of pH, as well as addressing common misconceptions regarding the use of indicators.

Equipment (per group)

- A4 print out of integrated instructions (laminated or within a plastic wallet)

All chemicals below to be supplied in 10 ml dropper bottles:

- Hydrochloric acid (0.1 mol dm⁻³), 12 drops
- Sodium hydroxide (0.1 mol dm⁻³), 12 drops
- pH 4 buffer solution, 12 drops
- pH 7 buffer solution, 12 drops
- pH 9 buffer solution, 12 drops
- Bromothymol blue solution (<1%), 20 drops
- Phenolphthalein (0.05%), 15 drops
- Methyl orange solution (<1%), 15 drops
- Universal indicator solution, 20 drops

Full preparation information and hazard classification is given below.

Safety equipment

- Eye protection: safety glasses to EN166 F

Method

Part 1 – testing indicators

Complete the experiment on the integrated instructions worksheet (slide 3).


Working down the columns:

1. Add 2 drops of strong acid to each square in the first column.
2. Add 2 drops of weak acid to each square in the second column.
3. Add 2 drops of neutral buffer to each square in the third column.
4. Add 2 drops of weak alkali to each square in the fourth column.
5. Add 2 drops of strong alkali to each square in the last column.

Working across the rows:

- Add 2 drops of bromothymol blue to each square in the top row.
- Add 2 drops of methyl orange to each square in the second row.
- Add 2 drops of phenolphthalein to each square in the third row.
- Add 2 drops of universal indicator to each square in the bottom row.

The expected results are shown in the image below.



Part 2 – making universal indicator

Complete the experiment on the integrated instructions worksheet (slide 4).

- Make the mixed indicator by adding 10 drops of bromothymol blue, 5 drops of methyl orange and 5 drops of phenolphthalein to a small beaker or vial.

Working down the columns:




- Add 2 drops of strong acid to each square in the column.
- Add 2 drops of weak acid to each square in the column.
- Add 2 drops of neutral buffer to each square in the column.
- Add 2 drops of weak alkali to each square in the column.
- Add 2 drops of strong alkali to each square in the column.







Working across the rows:




- Add 2 drops of your mixed indicator to each square in the top row.
- Add 2 drops of universal indicator to each square in the bottom row.

Safety

- Read our standard health and safety guidance (rsc.li/3zyJLkx) and carry out a risk assessment before running any live practical.
- Hazard classification may vary depending on supplier.

Chemical supplied for the practical	Preparation
Hydrochloric acid, 0.1 mol dm ⁻³ HCl(aq) Not currently classified as hazardous.	Hydrochloric acid, concentrated (≥6.8 M) HCl(l) MW = 35.46 g mol ⁻¹  DANGER Corrosive to skin and eyes. Irritant if inhaled. Refer to CLEAPSS Hazcard 47A. Prepare solution following CLEAPSS recipe sheet RB043.
Sodium hydroxide, 0.1 mol dm ⁻³ NaOH(aq) Not currently classified as hazardous.	Sodium hydroxide, solid NaOH(s) MW = 40.00 g mol ⁻¹  DANGER Corrosive to skin and eyes. Refer to CLEAPSS Hazcard 91A. Prepare solution following CLEAPSS recipe sheet RB085.
Buffer solutions pH 4, 7 and 9 Not currently classified as hazardous.	Buffer solutions Buffer tablets of pH 4, 7 and 9 can be used or alternatively prepare the buffer solutions following CLEAPSS recipe card RB018. If following the recipe card: Sodium tetraborate-10-water (borax)  DANGER May damage fertility and the unborn child. Refer to CLEAPSS Hazcards 13B, 09A and 14A.

<p>Phenolphthalein 0.05% solution</p>  <p>WARNING Flammable CLEAPSS Hazcard 032.</p>	<p>Phenolphthalein, solid DANGER</p>  <p>May include a risk of one or more hazard statements including, but not limited to: May be fatal/toxic if inhaled or in contact with skin or eyes. May cause breathing difficulties if inhaled. May damage fertility or the unborn child.</p> <p>Check CLEAPSS Hazcard 032 for most up to date guidance.</p> <p>Ethanol (and IDA) WARNING Flammable</p>  <p>CLEAPSS Hazcard 040A.</p> <p>Prepare solution following CLEAPSS recipe sheet RB046.</p>
<p>Methyl orange (<1%) solution</p>  <p>WARNING Flammable CLEAPSS Hazcard 032</p>	<p>Methyl orange, solid DANGER</p>  <p>May include a risk of one or more hazard statements including, but not limited to: May be fatal/toxic if inhaled or in contact with skin or eyes May cause breathing difficulties if inhaled May damage fertility or the unborn child</p> <p>Check CLEAPSS Hazcard 032 for most up to date guidance.</p> <p>Ethanol (and IDA) WARNING Flammable</p>  <p>CLEAPSS Hazcard 040A.</p> <p>Prepare solution following CLEAPSS recipe sheet RB046.</p>

<p>Bromothymol blue (<1%) solution</p>  <p>WARNING Flammable CLEAPSS Hazcard 032.</p>	<p>Bromothymol blue, solid Not currently classified as hazardous.</p> <p>Ethanol (and IDA) WARNING Flammable</p>  <p>CLEAPSS Hazcard 040A. Prepare solution following CLEAPSS recipe card RB046.</p>
<p>Universal indicator solution</p>  <p>WARNING Flammable CLEAPSS Hazcard 032.</p>	<p>Universal indicator solution</p> <p>Either use a commercial universal indicator or prepare a solution following CLEAPSS recipe card RB047.</p> <p>Follow the recipe in conjunction with the advice on CLEAPSS Hazcards 032 and 040A.</p>

Disposal

- Wipe down the laminated sheet/plastic wallet with a paper towel and dispose of the towel in laboratory waste.
- Rinse and dry the laminated sheet/plastic wallet.