



A microscale acid–base titration

Learning objectives

- 1 Safely carry out a microscale titration of sodium hydroxide and hydrochloric acid.
- 2 Use practical results to calculate an unknown concentration.

Introduction

In this experiment, you will use **microscale titration apparatus** to carry out an accurate titration on a much smaller scale. You will fill the microscale burette with a known concentration of hydrochloric acid and use a microscale pipette to transfer sodium hydroxide solution to a beaker. You will then carry out the titration, using your results to calculate the concentration of sodium hydroxide.

This technique can be a little fiddly at first! Microscale techniques allow us to work quicker and more safely (as less chemicals are used) while still maintaining accuracy in our results.

Equipment

Apparatus

- Graduated glass pipette, 2 cm³
- Pipette, 1 cm³, and pipette filler to fit (or a 1 cm³ plastic syringe)
- Plastic syringe, 10 cm³
- Fine-tip poly(ethene) dropping pipette
- Small lengths of rubber, plastic or silicone tubing
- Beakers, 10 cm³, x 2
- Clamp stand with two bosses and clamps
- Safety glasses


Chemicals

- Dilute hydrochloric acid, 0.10 M, about 10 cm³
- Sodium hydroxide solution, approx. 0.1 M (IRRITANT), about 10 cm³
- Phenolphthalein indicator solution (HIGHLY FLAMMABLE), a few drops





Method



1. Clamp the microscale burette as shown in the diagram. To fill the microscale burette, first push the syringe plunger completely down to ensure there is no air present inside. Place the tip of the microscale burette into the 0.10 M hydrochloric acid and slowly raise the plunger, making sure no air bubbles are drawn in. Fill all the way to the zero mark.
2. Use the 1 cm³ microscale pipette and pipette filler to transfer exactly 1.0 cm³ of the sodium hydroxide solution into a clean 10 cm³ beaker.
3. Add one small drop (no more!) of phenolphthalein indicator solution to the sodium hydroxide solution.
4. Adjust the position of the microscale burette so that the tip is just below the surface of the sodium hydroxide and indicator solution in the beaker.
5. Titrate the acid solution into the alkali by pressing down on the syringe plunger **very gently**, swirling to allow each tiny addition to mix and react before adding more.
6. Continue until the colour of the indicator just turns from pink to permanently colourless.
7. Record the volume of hydrochloric acid added at that point.
8. Repeat the titration until you get reproducible measurements – concordant results within 0.1 cm³ of each other.



Results

Draw a results table to record your data for this practical.

Calculating the concentration of sodium hydroxide solution

1. Write a balanced symbol equation for the reaction between sodium hydroxide and hydrochloric acid, and the ionic equation for the reaction.

-
-
-

-
-
-

-
- This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.