

# Some reactions of sulfur dioxide

## Instructions

1. Cover worksheet with a clear plastic sheet.
2. Place the base of the petri dish directly over the circle below. Place the reaction vessel in the centre.
3. Place two drops of full-range indicator solution where shown.
4. At another corner of the triangle place two drops of ammonia solution. Place the lid on the petri dish and wait for the indicator drop to change colour.
5. Remove the lid from the petri dish and, using a piece of tissue, mop up the drop of ammonia.
6. At the two remaining corners of the triangle add the two other test solutions.
7. Add a small quantity of sodium sulphite powder to the reaction vessel followed by three drops of hydrochloric acid. Quickly replace the lid.
8. Record all your observations over the next 15 min.

## Question

1. What explanations can you give for your observations?

## Health & Safety

Students must wear eye protection.

Sulfur dioxide is toxic and is a particular problem for asthmatics. Only use a very small amount of sulfite and acid to keep the sulfur dioxide production to a minimum. A risk-assessment should include any individual sensitivities.

Hydrochloric acid  $1 \text{ mol dm}^{-3}$ , Potassium iodide  $0.2 \text{ mol dm}^{-3}$ , Potassium iodate(V)  $0.1 \text{ mol dm}^{-3}$  and Potassium manganate(VII)  $0.01 \text{ mol dm}^{-3}$  are all of low hazard.

Sulfuric acid  $1 \text{ mol dm}^{-3}$  is a skin/eye irritant.

Sodium sulfite powder. Is a skin, eye and respiratory irritant.

Depending on its formulation full range indicator can still be flammable when diluted 1:1 with water. Keep away from sources of ignition.

## Credits

© Royal Society of Chemistry

*Health & safety checked May 2018*

Page last updated August 2018



